

Measurement of forces



Physics

Mechanics

Forces, work, power & energy



Difficulty level

easy



Group size

2



Preparation time

10 minutes




Execution time

10 minutes



Teacher information

Application



Experimental setup

In order to cause a constant acceleration \vec{a} at a mass m , a certain force \vec{F} is required according to Newton's second law:

$$\vec{F} = m \cdot \vec{a}$$

In practice, the force is often inferred from a known linear relationship between the acting force and an easily measured quantity. Examples are the deformation of an elastic material. In this case we use the length change of coil springs in so-called spring force meters, although this experiment does not specifically deal with Hooke's law.

The acceleration here is the acceleration of gravity $g = 9.81m/s^2$.

Other teacher information (1/2)

PHYWE
excellence in science

Prior knowledge



It is only assumed that students have a basic understanding of the concept of mass. The students should also know that the earth's gravitational field causes a permanent acceleration of gravity towards the centre of the earth.

The SI unit of force F is a Newton N :

$$1 N = 1 kg \cdot \frac{m}{s^2}$$

Scientific Principle



The students are to determine the weight forces of different objects via the deflection of a dynamometer.

Other teacher information (2/2)

PHYWE
excellence in science

Learning objective



The students should become acquainted with the dependence of the reading of spring balances on their orientation and with their adjustment..

Tasks



For this purpose, a spring balance adjusted downwards in the position of use shall be used to determine the weight of the following objects:

1. Helical spring
2. Holding pin
3. Support rod

Safety Instructions

PHYWE
excellence in science

The general instructions for safe experimentation in science lessons to be applied to this experiment.

PHYWE
excellence in science

Student Information

Motivation

PHYWE
excellence in science



Astronaut on the moon

As you know, every body can be assigned a mass m . This mass is independent of whether you are on the earth or the moon. So why is it that a man can jump much higher and farther on the moon than on the earth?

This is due to the local gravity field. This ensures that we are permanently accelerated towards the ground, because the so-called weight force is acting on us. This weight force is lower on the moon for an equal mass.

The weight force of a body, which results from its mass, is a fundamental quantity of physics. In this experiment you will learn how to determine it.

Tasks

PHYWE
excellence in science



In this experiment you will become familiar with the position dependence and the adjustment of spring balances.


With the spring balance adjusted downwards in the position of use, you can then determine the weight of various objects.

Material

Position	Material	Item No.	Quantity
1	Support rod with hole, stainless steel, 10 cm	02036-01	1
2	Helical spring, 3 N/m	02220-00	1
3	Spring balance,transparent, 1 N	03065-02	1
4	Spring balance,transparent, 2 N	03065-03	1
5	Holding pin	03949-00	1
6	Balance pan, plastic	03951-00	1

Set-up

PHYWE
excellence in science



Assembly of the balance pan

Assemble the balance pan as shown in the picture.

Procedure (1/6)

PHYWE
excellence in science




Spring balance vertical

Hold the spring balance (2 N) vertically so that the zero on the scale is at the top.

Observe the position of the marker in relation to the scale.

Procedure (2/6)

PHYWE
excellence in science




Spring balance horizontal

Now hold the spring balance (2 N) horizontally.

Observe the position of the marker in relation to the scale.

Procedure (3/6)

PHYWE
excellence in science



Spring balance vertical
(upside down)

Hold the spring balance (2 N) vertically so that the zero on the scale is at the bottom.

Observe the position of the marker in relation to the scale.

Procedure (4/6)

PHYWE
excellence in science



Adjusting the scale

- Adjust the scale of the spring balance (upside down / zero point of the scale below) by loosening the counter screw at the head end and turning the hook until the indicator mark is exactly at zero. Then tighten the screw again.
- Now hold the spring balance vertically downwards (zero point up), then horizontally and read the display each time. Enter the measured values in Table 1 in the report.

Procedure (5/6)


PHYWE
excellence in science



Adjustment to zero

- Adjust the spring balance 2 N held vertically downwards to zero (zero point of the scale at the top).

Procedure (6/6)

PHYWE
excellence in science

Measuring the helical spring and the support rod

- Hang the balance pan on the hook and note the measured value.
- Place the helical spring, the holding pin and the support rod on the balance pan one after the other.
- Note all measured values in Table 2 in the report.
- Repeat this series of measurements with the 1N spring balance and note these values in Table 2.

PHYWE
excellence in science

Report

Table 1



Spring balance position	$F [N]$
upside down	
horizontal	
vertical	

Enter your measurement results in the table.

Table 2



Enter your measurement results in the table. Use the values in the table to calculate the weight of the 3 objects without balance pan and enter the results.

Spring balance	Weight force $F [N]$			
	with weighing pan		without weighing pan	
	2 N	1 N	2 N	1 N
Balance pan				
Helical spring				
Holding pin				
Support rod				

Task 1

As you have noticed, the display changes depending on the position of the spring balance. It's true:

- The measured value is always highest for the normal vertical position.
- In the vertical position upside down the measured value is always the smallest.
- The measured value is always smallest for the normal vertical position.
- The value of the horizontal position is always between those of the two vertical ones.
- In horizontal position the scale always shows zero.

✔ Check

Task 2

Reasons for the deviations in the 3 different positions of the spring balances are:

- Optical distortions caused by the change in position relative to the eye.
- The measuring accuracy of the respective spring balances.
- The dead weight of the suspension on the spring.
- There are no deviations.
- The directional dependence of the earth's gravitational field.

✔ Check

Task 3

On the 1 N or 2 N spring balance, a graduation mark means

Always a change of 0.02 N.


Always a change of 0.1 N.

A change of 0.01 N or 0.02 N.

A change of 0.1 N or 0.2 N.

Slide	Score/Total
Slide 20: Position-dependent changes	0/3
Slide 21: Justification of the amendments	0/2
Slide 22: force gauge graduations	0/1

Total amount  0/6

 Solutions

 Repeat

 Exporting text